



BODY MAKEUP REMOVER

PAINLESS REMOVAL OF PROS-AIDE

The Body Makeup Remover (BMR) from Mavidon gently removes **Pros-aide®**, **collodion** and **most other adhesives** from your skin or hair for painless removal.

BMR is a proprietary blend of mild solvents designed to penetrate tapes and soften adhesives to facilitate removal. Frequently adhesive removers leave residue on the skin. With its slow evaporation rate BMR will leave less (if any) residue of adhesive. To remove any residue, simply wash off completely with any mild soap and water.

The product may be used for all age groups

BMR May also be used to clean (non-LDPE plastic) devices by soaking for 10-15 minutes. Wash any residue off with soap and water.

Intended Use:

Collodion: Apply liberally to dried collodion for 1-2 minutes. After soaking each site the collodion will be loosening up. The oily feeling is a moisturizing oil which we included in our product to help dissolve collodion and protects the skin and hair from drying out. However, please take caution because the remover will melt LDPE plastics (hearing aids, glasses, car seats, phones....etc)

Pros-aide: Apply liberally to dried collodion for 1-2 minutes. After soaking each site the collodion will be loosening up. Wash any residue with soap and water.

Body Paint: Simply apply BMR on the paint and wipe down. Wash any residue with soap and water.

Use with care:

- May discolor or remove furniture paints.
- May melt low density plastics.
- May dissolve exam gloves.
- For External Use Only!
- May Cause Eye Irritation.
- Do not use on open wounds!**
- Keep Container tightly closed.
- Discontinue if skin irritation occurs.

* Pros-aide is a registered trademark off ADM Tronics Unlimited, Inc.

Market: Professional Makeup Artist and Consumer Healthcare

Composition: Proprietary blend of mild solvents and moisturizers

Shelf life: 3 years from production date.

Product number and Packaging:

Standard sizes are:

- AT201-4 4 oz Bottles
- AT201-P 16 oz Bottles



Mavidon
6625 White Dr
Riviera Beach, FL 33407 USA
www.mavidon.com
561-585-2227

