



YOUR SOURCE FOR ATTACHMENT
AND REMOVAL OF EEG ELECTRODES

MORE THAN COLLODION!

MAVIDON SILENT AIR PUMP **HD**

HIGHER OUTPUT

MD0029

THE QUIET ONE

Instructions:

- Place on flat surface
- Attached air hose in the quick release adaptor.
- Plug in a 110 Watt, - grounded outlet.
- Turn pump on
- Press foot pedal and pump will automatically turn on.
- Release foot pedal to turn off pump.
- Use electrode applicator, supplied, to help dry collodion.



Part Number MD0029

Specs:

- Volts 120
- 60 HZ
- 2.4 Amps
- CFM 1.86
- PSI: Auto Shut-off: 40 psi
- Noise Level :50 db/A
- 11.5 lbs
- CSA/US Label
- Dimension 7 x 6 x 10
- On/off Switch

1.86 CFM

Warranty

- 1 year warranty on motor from date of purchase.

Distributed by: Mavidon® Riviera Beach, FL 800-654-0385
www.mavidon.com

Made in Italy